

Technology of the GeneXpert[®] System



Agenda

Cepheid System Families

Technology of the GeneXpert[®] System

Technology of the Xpert[®] Cartridge

Connectivity

Power & Safety Requirements

Cartridge disposal



Learning Objectives

The General Objective of this Module Is to Give You an Understanding of the Technology of the GeneXpert[®] System

At the end of the training, you will be able to:

- Recognize and recall the different GeneXpert Systems
- Explain the technology of the GeneXpert and how the cartridge works
- Summarize the basics of GeneXpert connectivity
- Recognize the GeneXpert System power requirement needs
- Recall basic safety precautions
- Explain the general disposal requirements of test cartridges


Polymerase Chain Reaction (PCR) Evolution



Cepheid Systems Family



History of Ongoing System Innovation



The Package

GeneXpert® System*

- Thermal and optical module(s)
- Computer System and GeneXpert Dx Software
- Barcode scanner



Cartridge

- Self-Contained
- Disposable
- Assay Definition Protocol (ADF)



Recommended Accessories

- UPS
- Surge protector

Optional Accessories

- Batteries/Power generator
- Printer


* FDA cleared

The Technology

The GeneXpert[®] System



Proven Technology




GeneXpert[®] Technology

- Integrated and closed system
 - No direct contact between instrument and sample to eliminate carry-over
 - The sample is enclosed in the cartridge
 - Integrated ultrasonic horn for cell lysis (when applicable)
- Fluid transfer: Micro fluidics based reconstitution, and automated filling
 - Advanced micro fluidics technologies to enable complex sample preparation processing protocols
 - Software driven motors for valve movement and integral hydraulic drives
- Multiple integrated controls to validate every step
- Automated protocol, data reduction and results interpretation

Automated Xpert® Protocol



GeneXpert[®] Module



"For complete details on the GeneXpert System, please refer to the GeneXpert System Operators Manual"


Cartridge Bay

I-CORE® Slit

PCR amplification & detection

Ultrasonic Horn

Lyses the sample (if applicable)



Plunger Rod


Facilitates movement of sample and reagents into different chambers

Valve Drive

Rotates the cartridge valve body to allow access to the different cartridge chambers

The I-CORE[®] Module

Building Block of the GeneXpert[®] System



Fan
Heater

Rapid, precise temperature control

Inserts into the I-CORE[®]




Cartridge
Sample preparation

Circuitry

Passes optical information to the computer for analysis and display

Optical Blocks

Optical analysis, detect and quantify up to 10 different DNA targets simultaneously



The Technology

Xpert Test Cartridge



Xpert[®] Cartridge

- Self-contained cartridge
- Avoids cross-contamination

Processing Chambers

Holds the sample, reagents, processed sample and waste solutions

Valve Body

Rotates and allows fluid to move to different chambers and to the reaction tube




Cartridge Lid

PCR Reaction Tube

Enables rapid thermal cycling and optical excitation/detection


Cartridge Foot

Dedicated for the Plunger



- Some chambers are used for sample processing
- Some chambers are for reagents


Internal View of the Cartridge



Cartridge Bottom View

Rotary valve is pre-aligned to fit with the valve in the module

- Do not rotate the rotary valve
- Alignment required to lock the module door and start the process




Bottom View of the Rotary Valve

PCR Reaction Chamber

Fluid transfer in the PCR Reaction Tube

- Do not touch the PCR reaction tube
- Wear gloves whenever you prepare a cartridge



Connectivity



Laboratory Information System (LIS) Compatibility

Result availability

- Improve Turn Around Time (TAT)
- Improve efficiency and effectiveness

Eliminate manual data entry

- Reduce risk of data entry errors
- Optimize workflow and simplify steps

Improve patient response and care

- Instant access to actionable results


The image shows two informational slides from Cepheid titled "LEAN: From Beginning to End".

The left slide is titled "How do you automate the most automated system in the molecular marketplace? Interface it to your LIS." It features a circular diagram with four quadrants: "Improve efficiency and effectiveness", "Reduce risk of manual entry", "Optimize workflow and simplify steps", and "Improve patient response and care". Below the diagram, it lists benefits such as "Reduce risk of manual entry", "Optimize workflow and simplify steps", "Improve TAT", and "Improve patient response and care". It also includes a quote from Bryan Desimone, Applications System Analyst, Legacy Health System, stating: "We have seen many workflow improvements after integrating the Cepheid to our LIS including increased productivity as well as a significant decrease in time to get our current error-free testing results."

The right slide is titled "Cepheid LIS Module" and includes a diagram showing the integration of the LIS, Cepheid, and the LIS. It lists features: "In-Direction Communication", "HL7 and/or ACTN compatible", "Automatic Data Transfer Capability", "Integration to Verification of Data", "Integration with Multiple Cepheid Systems at one time", and "Complete Interoperability for all Cepheid Systems". It also includes a quote from Bryan Desimone: "We've seen a great reduction in our turn-around times now that the built-in automatic interface to the instrument. Also helping reduce these times is that the results automatically get sent into our LIS which means we no longer have to type and fax either the results and then fax another 'back-entry sheet.'"

Both slides include a "Frequently Asked Questions" section with questions and answers regarding LIS compatibility, integration, and support.

Cepheid C360 Connectivity*



Power and Safety Requirements



System Power Consumption Requirements

- Power Supply: Auto Rating
- GeneXpert Energy Consumption information

System Size	On Power Mode Consumption (W)	Annual Energy Consumption (KWh)	Standby Power Consumption
GX-1	61	263	58
GX-II	85	372	71
GX-IV	100	489	83
GX-XVI	270	1168	170
Infinity-48	2426	5840	1248
Infinity-80	2426	5840	1248

- Computer Energy Consumption information
 - Laptop 350 kW
 - Desktop 350 kW

Safety Precautions

Injury may result if the instrument is not lifted properly.



It is possible for you or the system to be exposed to biological hazards.



The GeneXpert instrument's enclosure is designed to protect operators from electrical shock hazards.



Disposal

- Biological specimens, transfer devices, and used cartridges should be considered capable of transmitting infectious agents requiring standard precautions
- Follow your institution's environmental waste procedures for proper disposal
- These materials may exhibit characteristics of chemical hazardous waste requiring specific national or regional disposal procedures
- If national or regional regulations do not provide clear direction on proper disposal, the biological specimens and used cartridges should be disposed of per WHO [World Health Organization] medical waste handling and disposal guidelines



Technical Assistance

- Before contacting Cepheid Technical Support, collect the following information:
 - Product name
 - Lot number
 - Serial number of the System
 - Error messages (if any)
 - Software version and, if applicable, Computer Service Tag number
- Log your complaint online using the following link <http://www.cepheid.com/us/support> :
Create a Support Case



Thank You.



www.Cepheid.com