



GeneXpert[®] Edge Software

Cepheid Training Centre



302-0567 Rev. B May 2022

© 2022 Cepheid

IVD. In vitro diagnostic medical device. May not be available in all countries.

Agenda Title

- 1 Intended Use
- 2 Starting the System
- 3 GeneXpert Software
- 4 Managing User
- 5 GeneXpert Edge Software
- 6 Q & A

Training Objectives

At the end of the training, users will be able to:

- Configure the GeneXpert[®] software
- Configure the GeneXpert[®] Edge software
- Prepare a cartridge and run the test
- Report and print the various software generated results
- Perform basic troubleshooting

Intended Use

- The GeneXpert[®] Edge system is an in vitro diagnostic device intended for use with Cepheid Xpert[®] test kits. The GeneXpert[®] Edge system automates and integrates sample preparation, nucleic acid amplification, and detection of the target sequence in simple or complex samples using real-time Polymerase Chain Reaction (PCR). The system is designed for hands-off processing of patient samples (specimens) and provides both summarized and detailed test results data in tabular and graphic formats.
- The GeneXpert[®] Edge system is intended to be used by **laboratory professionals** or **specifically-trained healthcare users** in a **laboratory and near patient test settings** as specified in the Cepheid[®] Xpert[®] test instructions for use.

Software Installation

- The GeneXpert® Edge System comes with the GeneXpert® Software & GeneXpert® Edge software pre-installed.
- The GeneXpert® Edge software is a Simplified UI(User Interface) to run tests.
- It uses GeneXpert® DX Database.
- GeneXpert® Edge software 1.0 is exclusively for running and view tests.
- Key features included is a Video guidance for cartridges Preparation





Starting up the system

Starting the Software and Instrument

1. Put on a clean pair of gloves.

2. Turn on the GeneXpert® Edge Instrument.




3. Turn on the tablet computer. The Windows® screen will display.



4. Swipe up on the lock screen.

5. The Windows login screen will display. Touch the Password field to display the keyboard. Enter the password.



Starting the Software and Instrument continued

6. Touch the arrow button at the right of the password entry area to log in.



7. The Welcome Screen appears after the software initializes. Touch the TOUCH HERE TO BEGIN button.



8. The VIEW PREVIOUS TESTS button will display. The RUN NEW TEST button will appear on the Home screen within 3 minutes






GeneXpert[®] Software

GeneXpert[®] Software

- All configuration activities are performed on GeneXpert[®] Software prior to creating tests in GeneXpert[®] EDGE

- Create/Define Users
- Install ADF
- System configuration
- Access to Reports




➤ including Installation Qualification Report, System log


- LIS Settings

Important Reminders !

- **Create Test is disabled in GeneXpert Software**



- GeneXpert Software version is displayed at the bottom of the Check status screen






Managing Users

Edit user type configuration


Define access rights to the software

1



2

Edit user configuration by selecting or deselecting



Adding new users


Create an Administrator Account (“Admin”)

Create an Administrator Account (“Admin”)


Define the person in the lab that will have full access to the software

The first account to be created **must be an Admin account**. There can be more than one Administrator account per system

1 Setup → User Administration



2 Click Add




Create an Administrator Account (“Admin”)

3 Enter the Username and Password

4 Choose the User type (Admin)


5 Click **OK**



The screenshot shows the 'Add User' dialog box with the following fields and options:

- User Name:** admin1
- Full Name:** Administration User
- Password:** [masked with dots]
- Confirm Password:** [masked with dots]
- User Type:** Admin (selected in a dropdown menu)
- Buttons:** OK and Cancel (the OK button is circled in orange)

Orange arrows point from the numbered instructions to the corresponding fields: from instruction 3 to the User Name field, from instruction 4 to the User Type dropdown, and from instruction 5 to the OK button.




A smaller version of the 'Add User' dialog box is shown below the main one. It has the same fields as the main dialog. The 'User Type' dropdown is open, showing a list of options: Basic, Admin, and another option partially visible. A pink arrow points to the 'Admin' option in the list.




Create an Xpress Account

Create an Xpress Account

1 Setup → User Administration



2 Click Add



Create an Xpress Account

Password and username must be between 6-10 characters

3 Enter the Username and Password

The screenshot shows a dialog box titled "Edit User" with the following fields and options:

- User Name: xpress
- Full Name: Xpress User
- Password: [masked with dots]
- Confirm Password: [masked with dots]
- User Type dropdown menu:
 - Xpress (selected)
 - Basic
 - Detail
 - Admin
 - Xpress
- Buttons: OK (circled in orange), Cancel

4 Choose the User type (Xpress)


5 Click on OK

Editing User Profile

How to edit user profile

1

Select the user to be edited from the user administration list.




2

Click on **Edit**

3


Revise the information, as desired, and then click **OK** to save the changes and close the Edit User dialog box.



Removing Users


How to remove users

1 Click to select the user to be removed



2 Click on **Remove**

3 To remove additional users, repeat Step 1 and Step 2. If you are done removing users, click **Close**






Enabling Login for GeneXpert[®] Edge Users

Create an Xpress Account


Define the person (s) will only have access to the GeneXpert[®] Edge software. This user cannot access the GeneXpert software

1



2

Click on «Add»




Enable user log-in

The default setting does not require a login for GeneXpert[®] Edge users.


If login is required, then configure the system as such

1 Go to System Configuration > Xpress settings



2

Tick box to allow log-in for GeneXpert[®] Edge software



Host Communication Settings

Host Communication Settings

- To configure the GeneXpert[®] Edge software for a Laboratory Information System (LIS) interface connection, contact Cepheid Technical Support.
- See the Technical Assistance section in the Preface for contact information



GeneXpert[®] Edge 1.0 System Software

GeneXpert[®] Edge Software Overview


- **GeneXpert[®] Edge 1.0** software
- Requires **GeneXpert[®] Software** in the background
- Compatible with Windows 10
- Languages Supported: **English ONLY**
- GeneXpert[®] Edge and GeneXpert[®] folder shortcuts on Desktop
- Computer is Configured to Auto Launch GeneXpert[®] Edge Software
- Video Guidance for Cartridge Preparation included during create test process

Start-up

- Please wait a little – the Edge software will start automatically



If the software does not open automatically, double click the icon on the desktop







Create a test

Create a new test

1



2



Make sure that your module is available and connected

Create a Test

1

**RUN
NEW TEST**

2

Click YES if
Sample/Patient ID
barcode is available

Click NO and enter it
manually

3

Scan sample ID



Check information and confirm to
proceed

4

Please confirm if you have entered the correct Patient/Sample ID?
662698701197

Scan cartridge
barcode

5

Select the appropriate cartridge and press the trigger, as shown, to scan the barcode.



Create a Test

6

Please confirm that the selected Assay (Test) is correct?

Xpert MTB-RIF Assay G4

YES

NO

8



Load the cartridge and close the module door

7

Patient/Sample ID

3401551882096

Assay

Xpert MTB-RIF Assay G4



Please refer to the package insert for specific sample and cartridge preparation instructions

CONTINUE →

Watch the video and follow the cartridge preparation steps

Monitoring Progress

You can monitor the progress of your test. Time is displayed in minutes





View and print results

View results

Display shows: Patient/Sample ID, Assay and Start Date plus Time

1

**VIEW
PREVIOUS
TESTS**

2

Patient/Sample ID	Assay	Start Time
3401551882096	Xpert MTB-RIF Assay G4	24/10/18 10:49:48
3401551882091	DEMO_Xpert MTB-RIF Ultra	19/10/18 12:18:54
3401551882092	DEMO_Xpert MTB-RIF Ultra	15/10/18 16:49:50
3401551882094	DEMO_Xpert MTB-RIF Ultra	01/10/18 14:19:49
3401551882099	DEMO_Xpert MTB-RIF Ultra	18/09/18 06:22:11
3401551882021	DEMO_Xpert MTB-RIF Ultra	12/09/18 01:49:35
3401551882022	DEMO_Xpert MTB-RIF Ultra	05/09/18 00:12:38
3401551882023	DEMO_Xpert MTB-RIF Ultra	06/08/18 06:31:36
3401551882024	DEMO_Xpert MTB-RIF Ultra	06/08/18 06:25:10
3401551882025	DEMO_Xpert MTB-RIF Ultra	06/08/18 06:19:12

Navigation controls on the right: Up arrow, Down arrow, SELECT button, Down arrow, Down arrow.

Double-click on a result line to view a specific result

Print Results

3

Patient/Sample ID

12345a99

Cartridge S/N

7143473072

Assay

Xpert HIV Qual

Result

HIV-1 DETECTED

Start Time

06/04/18 09:05:23

Test Disclaimer

For In Vitro Diagnostic Use Only.

4

Ensure that your printer is connected and select print results

5

PRINT RESULT



Exiting the software


How to Exit Edge Software

Do NOT exit the software whilst a test is running

1



2



Select Yes to confirm and exit the software





Troubleshooting

NO RESULT/INVALID

Repeat the test by following the retest procedure

Patient/Sample ID	3401551882096	Cartridge S/N	235956652
Assay	Xpert MTB-RIF Assay G4		
Result	NO RESULT - REPEAT TEST		
		Start Time	24/10/18 10:49:48
		Test Disclaimer	For In Vitro Diagnostic Use Only.

[PRINT RESULT](#) 





Technical assistance

Technical Assistance

- Before contacting Cepheid Technical Support, collect the following information:
 - Product name
 - Lot number
 - Serial number of the System
 - Error messages (if any)
 - Software version and, if applicable, Computer Service Tag number
- Log your complaint using the following link <http://www.cepheid.com/us/support> - **Create a Support Case**

Region	Telephone	Technical Support Email
US	+ 1 888 838 3222	techsupport@cepheid.com
Australia and New Zealand	+ 1800 107 884 (AU) + 0800 001 028 (NZ)	techsupportANZ@cepheid.com
Brazil and Latin America	+ 55 11 3524 8373	latamsupport@cepheid.com
China	+ 86 021 5406 5387	techsupportchina@cepheid.com
France	+ 33 563 825 319	support@cepheideurope.com
Germany	+ 49 69 710 480 480	support@cepheideurope.com
India, Bangladesh, Bhutan, Nepal, and Sri Lanka	+ 91 11 48353010	techsupportindia@cepheid.com
Italy	+ 39 800 902 567	support@cepheideurope.com
South Africa	+ 27 861 22 76 35	support@cepheideurope.com
United Kingdom	+ 44 3303 332 533	support@cepheideurope.com
Other European, Middle East, and African countries	+ 33 563 825 319 + 971 4 253 3218	support@cepheideurope.com
Belgium, Netherlands and Luxembourg	+33 563 825 319	support@cepheideurope.com
Other European, Middle East, and African Countries	+ 33 563 825 319 + 971 4 253 3218	support@cepheideurope.com
Other countries not listed	+ 1 408 400 8495	techsupport@cepheid.com

Information on GeneXpert Edge System

- **GeneXpert[®] Edge User Manual**





Q&A



Thank You

www.cepheid.com