Xpert C.difficile

REF GXCDIFFICILE-CN-10

Cepheid AB
Röntgenvägen 5
SE-171 54 Solna
Sweden

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Cepheid AB
Röntgenvägen 5
SE-171 54 Solna
Sweden

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1. Larson HE, Price AB, Honour P, Borriello SP. *Clostridium difficile* and the aetiology of pseudomembranous colitis, Lancet 1978; 1:1063-1066.
2. Bartlett JG. Clinical practice. Antibiotic-associated diarrhea. N Engl J Med 2002; 31:334-339
3. Borriello SP. The influence of the normal flora on *Clostridium difficile* colonization of the gut. Ann Med 1990;22-61-7
4. Bignardi GE. Risk factors for *Clostridium difficile* infection. J Hosp Infect 1998; 40:1-15.
5. Kelly CP, Pothoulakis C, Lamont JT. *Clostridium difficile* colitis. N Engl J Med 1994; 330:257-262.
6. Braun VT, Hundsberger P, Leukel M, Sauerborn and C. von Eichel-Striber. Definition of the single integration site of the pathogenicity locus of *Clostridium difficile*. 1996; Gene 181:29-38.
7. Hammond GA, Johnson JL. The toxigenic element of *Clostridium difficile* strain VPI 10463. Microb Pathog. 1995;19:203-213.
8. Sambol SPM, Merrigan D, Lysterly DN, Gerding, Johnson S. Toxin gene analysis of a variant strain of *Clostridium difficile* that causes human clinical disease. Infect. Immun 2000;68:5480-5487.
9. Goncalves C, Decre D, Barbut F, Burghoffer B, Petit JC. Prevalence and characterization of a binary toxin (actin-specific ADP-ribosyl-transferase) from *Clostridium difficile*. J Clin Microbiol 2004;42:1933-9
10. Stubbs S, Rupnik M, Gibert M, Brazier J, Duerden B, Popoff M. Production of actin-specific ADP-ribosyltransferase (binary toxin) by strains of *Clostridium difficile*. FEMS Microbiol Lett 2000;186:307-12.
11. Popoff MR, Rubin EJ, Gill DM, Boquet P. Action-specific ADP-ribotransferase produced by a *Clostridium difficile* strain. Infect Immun 1998;56:2299-306.
12. Kuijper EJ, Coignard B, Tull P. ESCMID Study Group for *Clostridium difficile*; EU Member States; European Centre for Disease Prevention and Control. Emergence of *Clostridium difficile*-associated disease in North America and Europe. Clin Microbiol Infect. 2006; Oct;12 Suppl 6:2-18. Review.
13. Curry SR, Marsh JW, Muto CA, O'Leary MM, Pasculle AW, Harrison LH. *tdc* genotypes associated with severe *Tdc* truncation in an epidemic clone and other strains of *Clostridium difficile*, J Clin Microbiol. 2007 Jan;45(1):215-21. Epub 2006 Oct 11. Erratum in: J Clin Microbiol. 2007 Jun;45(6):2103.
14. MacCannell DR, Louie TJ, Gregson DB, Laverdiere M, Labbe AC, Laing F, Henwick S. Molecular analysis of *Clostridium difficile* PCR ribotype 027 isolates from Eastern and Western Canada, J Clin Microbiol. 2006 Jun;44(6):2147-52.
15. Wilkins TD, Lysterly DM. *Clostridium difficile* testing: after 20 years, still challenging. Clin Microbiol. 2003 Feb;41(2):531-4. Review.
16. Delmee M. Laboratory diagnosis of *Clostridium difficile* disease. Clin Microbiol Infect. 2001 Aug;7(8):411-6. Review.
17. Poutanen SM, Simor AE. *Clostridium difficile*-associated diarrhea in adults. CMAJ. 2004 Jul 6;171(1):51-8. Review.
18. Centers for Disease Control and Prevention. Biosafety in microbiological and biomedical laboratories. Richmond JY and McKinney RW (eds) (1993). HHS Publication number (CDC) 93-8395.
19. Clinical and Laboratory Standards Institute (formerly National Committee for Clinical Laboratory Standards). Protection of laboratory workers from occupationally acquired infections; Approved Guideline. Document M29 (refer to latest edition).

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Copyright © Cepheid. Cepheid®, Cepheid logo, GeneXpert®, 和 Xpert®均为 Ceph 商标。版权归 Cepheid 所有。