

# GeneXpert<sup>®</sup> Dx


## Reference Guide



## Good Laboratory Practice

Real Time Polymerase Chain Reaction (RT-PCR) is a standard laboratory testing method that is used to select a specific sequence of DNA or RNA. This segment is amplified exponentially which creates billions of detectable copies. This technique has become an important tool in clinical laboratories for the detection of infectious pathogens at extremely low levels. This highly sensitive technique makes RT-PCR highly susceptible to cross-contamination, particularly from sample to sample transfer, if proper clean molecular technique is not used. Implementation of safeguards and strict adherence to robust protocols is often sufficient to ensure that cross-contamination in the molecular laboratory is a rare event.

## Follow general CMS guidance for good laboratory practice



### [PDF] GOOD LABORATORY PRACTICES

<https://www.cms.gov/Regulations-and-Guidance/...>



**File Size:** 33KB

**Page Count:** 3

GOOD LABORATORY PRACTICES 1) Keep the manufacturer's product insert for the laboratory test in use and be sure it is available to the testing personnel. Use the manufacturer's product insert for the kit currently in ... +

# Preventing Cross Contamination

## Use of Personal Protective Equipment (PPE)

**Gloves:** Change gloves after touching a sample. The outside of the sample harbors much of the sample DNA/RNA that transfers to the surface of gloves.

**Lab coats:** Wear a lab coat while processing samples. Wearing lab coats will prevent the transfer of sample DNA/RNA to other areas of the room.

**Eye/face protection:** Wear surgical masks, face shields, or other physical barriers, like a splash shield for procedures with a high likelihood of generating droplets or aerosols.

## Cleaning

**Bleach:** Use a final concentration of 1:10 dilution of 5% household chlorine bleach (used within 1 day of preparation). Final active chlorine concentration should be 0.5%.

**70% ethanol:** Use only 70% ethanol or denatured ethanol (70% ethanol containing 5% methanol and 5% isopropanol).

- Disposable lint-free wipes
- Disposable paper towels

## Reagent Storage

Store reagents according to their expected storage conditions in the Information for Use. In addition, cartridges should be kept in their original boxes with the lid shut.

## Sample Setup

**Dirty area (work area):** Area where samples and controls are processed.

**Clean area (loading area):** Area where the prepared cartridge is loaded onto the instrument.

## Cartridge Disposal

Used cartridges may contain potentially infectious materials, as well as highly amplified PCR target(s). **Do not open or attempt to alter any part of the cartridge for disposal.**

Each state has different regulations for classifying regulated medical waste (RMW). The first step to safe biohazard waste disposal is to check with your state's Department of Health to learn the specific regulations you'll need to follow.

## Maintenance

Instrument maintenance is required to be performed according to the user guide or operator manual. Some of the maintenance is described in this reference guide, however not all requirements are covered.


## Starting up the system

See the GeneXpert Dx System Operator Manual or Assay Instructions for Use for detailed information..


1. Turn the power switch on the instrument to the **ON** position.  
The blue light on the front panel will light up.



2. Turn the computer ON.
3. User-Account: **Cepheid-Admin**  
Password: **cphd**



4. The GeneXpert Dx software starts automatically. Enter user name and password if applicable.
5. In the **Check Status** screen, verify that all the modules are **available**.




Modules						Tests Since Launch							
Module Name	Assay	Sample ID	Progress	Status	Remaining Test Time	Sample ID	Mod Name	User	Result	Assay	Status	Error Status	Start Date
A1			Available										
A2			Available										
A3			Available										
A4			Available										

# Shutting down the system

**Note: Restart the system once per week.**

**When performing this task, make sure no tests are running.**


1. Exit the GeneXpert Dx software.



2. Turn the computer OFF through the Windows home button.



3. Turn the power switch on the instrument to the **OFF** position.  
The blue light on the front panel will turn off.



**Note: Wait 2 minutes before restarting the system.**

## Common GeneXpert Dx Menus

See Appendix A of the Operator Manual for the complete list



### User

- Login
- Change Password
- Logout
- Exit

### Data Management

- Archive Test
- Retrieve Test

### Reports

- Specimen Report
- Patient Report
- Patient Trend Report
- Control Trend Report
- System Log
- Assay Statistics Report
- Installation Qualification

### Setup

- User Administration (Create/Edit Users)
- User Type Configuration
- System Configuration
- Assign Instrument Letter

### Maintenance

- Module Reporters
- Plunger Rod Maintenance
- Valve Maintenance
- Perform Self-Test
- Open Module Door or Update EEPROM
- Exclude Modules from Test command

### About


- About GeneXpert Dx System

# Creating A Test

1. Click on **Create Test** from the main menu of the GeneXpert® Dx




2. Enter or scan the Sample ID and Patient ID (if applicable).  
Scan the barcode on the cartridge.



3. Enter or verify the correct information for the following sections (if applicable):

- Patient ID
- Sample ID

A screenshot of the 'Create Test' dialog box in the software. The dialog has several input fields and dropdown menus. Two orange arrows point from the text in step 3 to the 'Patient ID' and 'Sample ID' fields. The 'Patient ID' field contains 'Patient ID' and the 'Sample ID' field contains '234567'. Other fields include 'Name' and 'Version' (dropdowns), 'Select Assay' (dropdown with 'Xpert Assay 1'), 'Select Module' (dropdown with 'A2'), 'Reagent Lot ID', 'Expiration Date' (YYYY/MM/DD), 'Cartridge S/N', 'Test Type' (dropdown with 'Specimen'), 'Sample Type' (dropdown with 'Other'), and 'Other Sample Type'. A 'Notes' text area is at the bottom. At the bottom of the dialog, there are three buttons: 'Start Test' (circled in orange), 'Scan Cartridge Barcode', and 'Cancel'.

4. Click on **Start Test** to begin test.

5. Load the cartridge in the module with the blinking green light.  
Close the module door until the green light stops blinking.






# Loading Assay Definition File (ADF)

**Note:** Importing of the Assay Definition File (ADF), located in the kit, is required only when adding a new assay for the first time or when an assay has been updated.

1. Insert the assay definition CD, located in the kit, into the computer's DVD drive.




2. Click **Define Assays**.




3. Click on **Import**.

4. Select the DVD drive.



5. Select **GeneXpert Systems** folder.

6. Select the .gxa file.




7. Click on **Import**.



# View Results and Generate/Print a Report

Click on **View Results**.





Click on **Report**.

Click on **View Test**.

Click on the checkbox(es) of the report to be printed or generated.

Double-click on the test to be viewed.



**Test Report**

Patient ID: 01000234390J7  
 Patient ID 2:  
 Sample ID: 15002496KR  
 Test Type: Specimen  
 Sample Type:

Assay Information

Assay	Assay Version	Assay Type
Xpert MTB-RIF Assay G4	5	In Vitro Diagnostic

Test Result: **MTB DETECTED LOW**  
**Rif Resistance NOT DETECTED**


Test and Analyte Result

Analyte Name	CI	EndPt	Analyte Result	Probe Check
Probe D	27.6	193	POS	PASS
Probe C	27.5	165	POS	PASS
Probe E	28.4	146	POS	PASS
Probe B	28.0	132	POS	PASS
SPC	28.0	277	NA	PASS
Probe A	26.8	135	POS	PASS
QC-1	0.0	0	NEG	PASS
QC-2	0.0	0	NEG	PASS

User: Mikobakteriologias Lab  
 Status: Done  
 Expiration Date: 15/05/16  
 S/W Version: 4.4a  
 Cartridge S/N: 333603671  
 Reagent Lot ID: 08202


Start Time: 26/01/15 11:17:08  
 End Time: 26/01/15 12:57:46  
 Instrument S/N: 702365  
 Module S/N: 638095  
 Module Name: A4

Notes:  
 Error Status: OK




# Patient ID Report (if applicable)

## 1. Select Patient Report



## 2. Enter the patient ID.



## 3. Click Preview PDF..



### Patient Report

Found Patient ID #2 = H112874895762R


- 2 Test(s) Found -

Patient ID: H112874895762R  
Sample ID: SD142231  
Assay: Xpert CDIFFICILE  
Assay Version: 3  
Test Result: **NEGATIVE**  
Start Time: 06/09/16 12:38:42  
Test Type: Specimen  
User: Detail User  
Status: Done  
Notes:

Patient ID: H112874895762R  
Sample ID: SD142231  
Assay: Xpert BCR-ABL Monitor IS  
Assay Version: 1  
Test Result: **ERROR**  
Start Time: 06/09/16 12:41:13  
Test Type: Specimen  
User: Detail User  
Status: Aborted  
Notes:



# Archiving and Purging

1. Select **Data Management** and **Archive Test**.




2. Highlight the tests to be archived. Click **Select Highlighted**, then click **OK**.

**Note:** Check **Purge** to remove archived tests from the database.




3. Click **Proceed**.




The file name is generated automatically.

4. Click **Save**.



5. Click **OK**.



The archived file can be found in the folder C:\GeneXpert\export

**Note:** If **Purge Selected Tests** was checked, confirm the selection by clicking **Yes**.




6. Copy archived data file to an external location.

# Cartridge Bay and Plunger Rod Cleaning

## Required Materials

- 1:10 dilution of household chlorine bleach prepared within the same day.  
Final Active Chlorine concentration should be 0.5%, regardless of the household bleach concentration in your country
- 70% ethanol or denatured ethanol (70% ethanol containing 5% methanol and 5% isopropanol)
- Lint-free wipes

1. Remove cartridge(s) from the module(s).
2. Click on **Maintenance** on the Menu Bar, select **Plunger Rod Maintenance**.
3. Select the module(s) to be cleaned and then select **Clean** or **Clean All**.
4. Click **OK**.



5. The plunger rod(s) in the selected module(s) lower(s) into the cartridge bay(s).
6. To clean:
  - A. Thoroughly moisten a lint-free wipe with a 1:10 solution of household chlorine bleach.
  - B. Vigorously wipe the plunger rod with the lint-free wipe. Using the same lint-free wipe, wipe the walls, ceiling, corners and edges of the cartridge bay, then wipe the inside of the door and the top lip of the door and discard the lint-free wipe.
  - C. Wait 2 minutes after wiping with the bleach solution.
  - D. Repeat steps A-C twice more, using a new lint-free wipe each time.
  - E. Wait 2 minutes after wiping with the bleach solution.
  - F. Thoroughly moisten a lint-free wipe with the 70% ethanol solution.
  - G. Repeat step B.
7. Once cleaning is completed, click **Move Up**.
8. Click **Close**.



Refer to the Operator Manual for additional Maintenance requirements/tasks.





**m**power

The molecular revolution is here.

